

SURGICAL PROCEDURE PACKS



01.

V.kit

Designed **with** Vets, **for** Vets

Procedure Packs

Pennine work in partnership with individual practices/surgeons to create customer specific, surgical theatre packs for a wide range of specialties.

We work closely with veterinarians to design their custom kit's which can contain a vast range of components such as; swabs, drapes, sutures, blades, tubing, suction devices, wound dressings, bandages, gloves, and gowns.

These are then arranged into the sterile procedure pack, ready for use.

These packs are hugely beneficial for a variety of reasons;

- Reduces lay-up times for procedures
- Easier to calculate the on-cost to the customer
- Reduces the potential for infection
- Reduces ordering and administration costs
- Ease of use
- Reduces emergency surgery response times
- Eradication of errors in setting up

Why Choose Pennine?

Having been one of the instigators of the production of procedure packs within the UK in 1992, it has since become a substantial part of Pennine's business and our competencies in this area are second to none.


Pennine Veterinary are one of the only pack manufacturers who offer truly bespoke procedure packs to the veterinary market. We enable our customers to choose every pack component in detail, even down to the drape material!

A Key Account Manager will have a portfolio of products and will walk you through the process, making as many alterations as you need along the way.

We provide your choice of sterile or non sterile samples to make sure you are totally happy with the finished pack.

Our flexible ordering options (standing orders, call off orders, or make-to-orders) allow stock holding to be tailored to the practice needs.

Procedure Packs



BRO-11_Rev 1_04/23

02.

Create Your Own

Step 1
Select your
components

We have a vast portfolio of items to choose from, allowing us to offer truly bespoke packs.

Step 2
Select your
drapes & gowns

We understand the surgeon's individual preferences for draping, and therefore provide a wide range of drape materials. We have additional options with fenestrations, adhesives, and reinforcements.

Step 3
Trial your pack
for free


Once we receive your desired specification, we will send a no obligation quote and sterile sample.

Any amendments can be made after the trial and further quotes will be sent.

Step 4
Sign off
& order

The lead time will vary depending on availability of the components chosen, but we will keep the customer updated along the way. From then on, we recommend using our standing order service to ensure packs will automatically be delivered as required, without the need to remember to re-order.

Pennine Veterinary also offer a variety of Standard Packs.



Surgical Packs

Our pre-made surgical packs are a great option for the practices that don't require large numbers of packs, but want a quick, 'off the shelf' solution to surgical consumables.

These are preassembled, sterile, and readily available to order.

Practices that have not used packs before may also benefit from these as a good introduction to procedure kits. They can act as a base to customise your own pack.

Orthopaedic Pack

Product Code	Description
PNV-00100	Orthopaedic Pack
1	Trolley cover 152cm x 190cm
1	Bowl 1000ml
1	Connecting tube, soft flexi rib 3m
4	Drape 75cm x 75cm with adhesive edge
1	Extremity drape 100cm x 160cm with 4cm elasticated fenestration
1	Kidney dish 800ml
5	Swabs 7.5cm x 10cm X-RAY
5	Lap swabs 30cm x 30cm X-RAY
1	Bulb Syringe

General Surgery Pack

Product Code	Description
PNV-00400	General Surgery Pack
1	Trolley cover 152cm x 190cm
1	Bowl 1000ml
1	Connecting tube, soft flexi rib 3m
4	Drape 75cm x 75cm with adhesive edge
1	Drape 150cm x 180cm with adhesive edge
1	Kidney dish 800ml
5	Swabs 7.5cm x 10cm X-RAY
5	Lap swabs 30cm x 30cm X-RAY
1	Bulb Syringe



We're here to help, please get in touch via telephone on +44 (0)1332 794880 or email us at customer.services@penninehealthcare.co.uk

04.

Surgical Packs

Soft Tissue Pack

Product Code	Description
VKT-09100	Soft Tissue Pack
1	Trolley cover 150cm x 190cm
4	Drape 75cm x 75cm with adhesive edge
1	Syringe 10ml luer slip
1	Kidney dish 800ml
2	Swab 10cm x 7.5cm X-RAY
1	Co-flex Bandage 5cm x 4.5cm
1	Blade size 11

Neuro Pack

Product Code	Description
VKT-09200	Neuro Pack
1	Trolley cover 150cm x 190cm
2	Drape 90cm x 75cm with adhesive edge
1	Drape 260cm x 150cm with adhesive edge
1	Drape 180cm x 175cm with adhesive edge
2	Drape tape 10cm x 50cm
1	Diathermy Bipolar forceps, with Lead
1	Connecting tube, 3m
1	Syringe 20ml luer slip
2	Swab 10cm x 10cm
1	Blade size 10
1	10fr Frazier sucker
1	Kidney dish 800ml



05.

Drape Packs

Our standard drape packs are extremely useful as an easy and good practice draping solution for varied animal sizes.

The adhesive edges allow you to create the necessary size of fenestration, without any need to cut into drapes.

There are four drapes in each pack to achieve the 'four square' drape system. The trolley cover wrapping creates a sterile field so the pack is ready to be used.

These drape packs are great value for money and save on ordering and set up times.

Large Drape Pack

Product Code	Description
VKT-04601	Large Drape Pack
1	Trolley cover 150cm x 190cm, zone reinforced, absorbent
4	Drape 90cm x 90cm with adhesive edge
1	Hand towel 38cm x 40cm

Small Drape Pack

Product Code	Description
VKT-01400	Small Drape Pack
1	Trolley cover, 150cm x 190cm, zone reinforced, absorbent
4	Drape 90cm x 75cm, Bi-laminate with adhesive edge



We're here to help, please get in touch via telephone on +44 (0)1332 794880 or email us at customer.services@penninehealthcare.co.uk

Spay Packs

Our Spay Packs are an excellent solution for spay surgeries as they are great value for money and have all the necessary consumable components in one sterile pouch, ready to use.

One kit means less time opening individual items, less waste, and a clear cost.

These are a stocked item, readily available to order.

Pennine Veterinary can also offer fully customised packs. Please see our pack brochure for further details.

Cat Spay Drape Pack

Product Code	Description
VKT-04400	Cat Spay Drape Pack
1	Blue wrap 60cm x 60cm
1	Drape 45cm x 50cm with 3cm x 5cm adhesive fenestration
2	Hand towels 38cm x 40cm
1	Scalpel size 15
5	Swab 10cm x 10cm X-RAY

Bitch Spay Drape Pack

Product Code	Description
VKT-04500	Bitch Spay Drape Pack
1	Blue wrap 60cm x 60cm
1	Drape 100cm x 100cm with 12cm x 9cm adhesive fenestration
2	Hand towels 38cm x 40cm
1	Scalpel size 15
10	Swab 10cm x 10cm X-RAY

Laparoscopic Bitch Spay Pack

Product Code	Description
VKT-07201	Laparoscopic Bitch Spay Kit
1	Drape, patient 300cm x 160cm, bi-laminate, non-reinforced, non-adhesive
1	Trolley cover, 150cm x 190cm, zone reinforced, absorbent
1	Insufflation set, with push fit / male luer connection, 3m x 6mm, square filter
1	Blade, size 11
1	Syringe, 20ml, luer slip
1	Needle, hypodermic, 21g, 1.5 inch, green
1	OR Towel, blue, x-ray, absorbent, 45 x 60cm
1	Bag, breather, 35cm x 50cm, 2fs
1	Suture, biosyn, 3-0, violet, 75cm, 1/2 circle round bodied needle
1	Bag, extras wallet, 18cm x 24cm
10	Swab, gauze, x-ray, 10cm x 10cm, 16ply, without tapes, white, 17 thread





Pennine Veterinary
Ivor Shaw Limited t/a Pennine Healthcare
300 City Gate, London Road, Derby, DE24 8WY, UK.
Tel: +44 (0)1332 794877
www.pennineveterinary.co.uk