

UROLOGY





ABOUT PENNINE


Pennine Healthcare is proud to be one of the UK's leading, employee owned, manufacturer and distributor, of single use, sterile Ward & Theatre Consumables and Custom Procedure Packs, as well as being one of the leading suppliers of PPE into healthcare settings. We market our products and services in over 50 countries worldwide and employ around 200 people at our purpose-built premises in Derby.

ABOUT OUR CATHETERS

Here at Pennine Healthcare, we understand the need for a reliable product range, to fulfil the need for those who need to empty their bladder via a catheter. Whatever your requirement, our catheters come in a variety of sizes and lengths to fulfil your individual needs.



NELATON (MALE) CATHETERS



Pennine's Intermittent Urethral Catheters are a single use, flexible tube inserted through the urethra into the bladder for the intermittent/transient passive evacuation of urine, which is then immediately removed.

- 39cm length (with connector)
- Medical grade PVC
- Sterile, individual peel pouches
- 100 units per inner polythene bag
- Supplied single wrapped (SW) as standard
- Latex & DEHP Free

Single Wrapped Product Code	NHS NPC Code	Qty Per Box	CH/Diam
NC-1206/SW	-	500	6ch ø 2.00mm
NC-1206/SW/1	-	100	6ch ø 2.00mm
NC-1208/SW	FST1000	500	8ch ø 2.66mm
NC-1208/SW/1	FKC180	100	8ch ø 2.66mm
NC-1210/SW	FST1002	500	10ch ø 3.33mm
NC-1210/SW/1	FKC184	100	10ch ø 3.33mm
NC-1212/SW	FST1004	500	12ch ø 4.00mm
NC-1212/SW/1	FKC186	100	12ch ø 4.00mm
NC-1214/SW	FST1006	500	14ch ø 4.66mm
NC-1214/SW/1	FKC194	100	14ch ø 4.66mm
NC-1216/SW	FST1007	500	16ch ø 5.33mm
NC-1216/SW/1	FKC196	100	16ch ø 5.33mm
NC-1218/SW	FST1008	500	18ch ø 6.00mm
NC-1218/SW/1	FST1009	100	18ch ø 6.00mm

Double Wrapped Product Code	NHS NPC Code	Qty Per Box	CH/Diam
NC-1210/DW	FKC183	400	10ch ø 3.33mm
NC-1212/DW	FKC185	400	12ch ø 4.00mm
NC-1214/DW	FKC193	400	14ch ø 4.66mm
NC-1216/DW	FKC197	400	16ch ø 5.33mm



For more information on our Catheter range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk

Please note: All dimensions are nominal.
Please note any products featuring a "*" indicates the product is distributed by Pennine.
For Regulatory and Technical information, Technical Data Sheets, are available upon request.


FEMALE CATHETERS





Pennine's Intermittent Urethral Catheters are a single use, flexible tube inserted through the urethra into the bladder for the intermittent/transient passive evacuation of urine, which is then immediately removed.

- 23cm length
- Medical grade PVC
- Sterile, individual peel pouches
- 100 units per inner polythene bag
- Latex & DEHP Free

Single Wrapped Product Code	NHS NPC Code	Qty Per Box	CH/Diam
FC-1406	-	1000	6ch ø 2.00mm
FC-1406/1	-	100	6ch ø 2.00mm
FC-1406/18/1	-	100	6ch ø 2.00mm
FC-1408	FST968	1000	8ch ø 2.66mm
FC-1408/1	FKC171	100	8ch ø 2.66mm
FC-1408/18/1	-	100	8ch ø 2.66mm
FC-1410	FST969	1000	10ch ø 3.33mm
FC-1410/1	FKC173	100	10ch ø 3.33mm
FC-1410/18/1	-	100	10ch ø 3.33mm
FC-1412	FST970	1000	12ch ø 4.00mm
FC-1412/1	-	100	12ch ø 4.00mm
FC-1412/18/1	-	100	12ch ø 4.00mm
FC-1414	FST975	1000	14ch ø 4.66mm
FC-1414/1	FST980	100	14ch ø 4.66mm
FC-1414/18/1	-	100	14ch ø 4.66mm
FC-1416	FST985	1000	16ch ø 5.33mm
FC-1416/1	FKC175	100	16ch ø 5.33mm
FC-1416/18/1	-	100	16ch ø 5.33mm
FC-1418	FST990	1000	18ch ø 6.00mm
FC-1418/1	-	100	18ch ø 6.00mm
FC-1418/18/1	-	100	18ch ø 6.00mm



URETHROTECH UCD®



The NICE-approved Urethral Catheterisation Device (UCD®) provides a safe second-line catheterisation solution when standard Foley-catheterisation has failed.

The integrated non-traumatic UCD-guidewire navigates the UCD-catheter safely into the bladder, preventing urethral injury particularly in men with an enlarged prostate. With the UCD® multiple traumatic unsuccessful catheterisation attempts are avoided, providing safer and cost-effective patient care.

- Safe and easy for patient and provider
- Ready-to-use consumable
- Suitable to be used by qualified doctors and nurses
- Improved patient care pathway
- Avoid invasive specialist interventions and associated costs
- Complications are reduced, providing safer and better care
- Clinically proven and evaluated by NICE*
- Purpose-designed and regulated medical device
- FDA 510(k), CE-mark, ISO 13485:2016
- Integrated hydrophilic Nitinol guidewire
- 100% Silicone 3-way Foley Catheter 16F
- The UCD-guidewire is removed after successful catheter placement
- Suitable for bladder irrigation if required
- Attached plug to close off guidewire (irrigation) channel

* Available for the UK market only



We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk



For more information on our Catheter range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk



Ivor Shaw Limited t/a Pennine Healthcare
300 City Gate, London Road, Derby, DE24 8WY, UK.
Tel: +44 (0)1332 794880
www.penninehealthcare.co.uk

 @arma_health  Pennine Healthcare  @Pennine_Health

BRO-6_Rev_4 03/2023