

SURGICAL CONSUMABLES



01.


Laparoscopic Surgery

A range of products specifically designed for use in laparoscopic surgical procedures.

Veress Needle

The Veress Needle is designed to introduce or remove gas from the peritoneal cavity as a therapeutic or surgical/radiological procedural method.

Product Code	Description
AELAP0001	120mm Veress Needle
AELAP0010	150mm Veress Needle



Camera Drape

The disposable camera sleeve is a sterile single-use drape made of tough polyethylene. The drape is the quickest and simplest way to create a sterile field for use with video/camera systems during operative laparoscopic procedure.

Product Code	Description
E6222	Camera cover



Insufflation Tubing Sets

Insufflation Tubing Set with Gauze Swab.

Product Code	Description
PEN-1250	Male Luer and Wolf Connector
PEN-1330E	Male Luer and Funnel Connector
PEN-1330LL	Male Luer and Male Luer
PEN-1331	Male Luer and 15/22mm Filter Adaptor
PEN-1334	Female Luer and Male Luer




We're here to help, please get in touch via telephone on +44 (0)1332 794880 or email us at customer.services@penninehealthcare.co.uk

Surgical Consumables

Cautery / Diathermy Pencil and Tip Cleaner

This Diathermy Pencil is a disposable electrosurgical pencil. Detachable tip electrodes, switch and buttons designed for smooth activation and supplied in sterile peel pouches.




Product Code	Description
PS555	Diathermy Pencil
TP4711/S	Disposable Adson Bipolar forceps with flying lead 15mm and 0.5 tip

Irrigation Bulb Syringe

Irrigation Bulb Syringe provides control of drawn or expelled volume by way of the design of its soft pliable bulb.

Ideal for use in clinical settings where a sterile element is not necessary, this sterile disposable irrigation syringe facilitates easy gripping. Designed to stand on its end, this medical bulb syringe helps prevent leakage to help you maintain a clean environment. This product has a volume of 60ml.



Product Code	Description	Box Qty
DYNDE20125	Bulb Irrigation Syringe 60ml with protective tip cover	50



We're here to help, please get in touch via telephone on +44 (0)1332 794880 or email us at customer.services@penninehealthcare.co.uk



03.

Frazier Suction Probes

These Frazier Probes are sterile single-use probes. With 'VC' vacuum control hole for precise control over the amount of suctioning. Used in open surgical procedures where delicate suctioning is required.


Product Code	Length (cm)	Size (ch)
AEMF1108	23	8
AEMF1110	23	10
AEMF1112	23	12



Micro-Flex Suction Probes

Our 'Micro-Flex Suction Probes are made from PVC and have a fingertip vacuum control. Designed for fine suctioning and in applications where controlled suction is required.

Product Code	Description
YS-3108	Microflex Suction VC 10fg 22.3cm
YS-3128	Microflex Suction VC 12fg 21.3cm
YS-3088	Microflex Suction VC 8fg 22.3cm



Feeding Tubes

A feeding tube is a small, soft, plastic tube placed through the nose or mouth into the stomach. Closed distal end with two lateral eyes to reduce risk of blockage.

Product Code	Length (cm)	CH/Tip
AEFT190540	40	5ch ø 1.7mm
AEFT190640	40	6ch ø 2.0mm
AEFT190840	40	8ch ø 2.7mm
AEFT191040	40	10ch ø 3.3mm



Esophageal Stethoscope

Esophageal Stethoscopes are used to ascertain the condition of the heart and lungs by listening to their action.

We have short (55cm) and long (110cm) sizes available.

55cm Esophageal Stethoscope

Product Code	Description
ESV-5024	55cm Esophageal Stethoscope fg24
ESV-5018	55cm Esophageal Stethoscope fg18
ESV-5012	55cm Esophageal Stethoscope fg12



110cm Esophageal Stethoscope

Product Code	Description
ESV-1124	110cm Esophageal Stethoscope fg24
ESV-1118	110cm Esophageal Stethoscope fg18
ESV-1112	110cm Esophageal Stethoscope fg12




We're here to help, please get in touch via telephone on +44 (0)1332 794880 or email us at customer.services@penninehealthcare.co.uk



Suction Connecting Tubing

Our Suction Connecting Tubing is a complete system to connect suction sources to suction waste collection systems, suction catheters, yankauers, suction probes and other suction devices.

Product Code	Description
CT-4012	Suction Connecting Tube 6mm x 1.8m F/F
CT-4022	Suction Connecting Tube 6mm x 2m F/F
CT-4024	Suction Connecting Tube 6mm x 2m F/F/M
CT-4034	Suction Connecting Tube 6mm x 3m F/F/M
CT-4034/VC	Suction Connecting Tube 6mm x 3m F/F/M VC
CT-4042	Suction Connecting Tube 6mm x 2.7m F/F
CT-4732	Suction Connecting Tube Flexi Rib 7mm x 3m F/F



Capped Vacuum Control Connector



The vacuum control unit is comprised of two PVC male connectors with a hole in between that can be covered by a finger during use to control the vacuum and suction flow rate of a device while in use. The vacuum control unit also has a cap that can be placed to cover the hole to allow for continuous suction without the need of the user to cover the hole with their finger.

Product Code	Type	Box Qty
VCU-2	Sterile	100




We're here to help, please get in touch via telephone on +44 (0)1332 794880 or email us at customer.services@penninehealthcare.co.uk

Link Yankauers & Sets

Our Link Yankauers are rigid, transparent suction devices for use with suction connecting tubing. For general purpose and fine suctioning, as well as suctioning of large volumes. Also available as a set, pre-attached to 2 or 3 metre tubing.

Product Code	Description
LY-3600	Yankauer Link Midi Tip 18fg 240mm
LY-3601	Yankauer Link Midi Tip VC 18fg 240mm
LY-3700	Yankauer Link Maxi Tip 22fg 275mm
LY-3701	Yankauer Link Maxi Tip VC 22fg 275mm
LY-3800	Yankauer Link Maxi Plus Tip 22fg 280mm
LYS-5620	Yankauer Link Set, Tube 6mm x 2m, Yankauer Midi LY-3600
LYS-5721	Yankauer Link Set, Tube 6mm x 2m, Yankauer Maxi VC LY-3701



We're here to help, please get in touch via telephone on +44 (0)1332 794880 or email us at customer.services@penninehealthcare.co.uk





Pennine Veterinary
Ivor Shaw Limited t/a Pennine Healthcare
300 City Gate, London Road, Derby, DE24 8WY, UK.
Tel: +44 (0)1332 794877
www.pennineveterinary.co.uk