

SURGICAL PROCEDURE PACKS



01.

V.kit

Designed **with** Vets, **for** Vets

About Pennine Veterinary

Having been one of the instigators of the production of procedure packs within the UK in 1992, it has since become a substantial part of Pennine's business and our competencies in this area are second to none.

Pennine Veterinary are one of the only pack manufacturers who offer truly bespoke procedure packs to the veterinary market. We enable our customers to choose every pack component in detail, even down to the drape material!

Our Key Account Manager will have a portfolio of products and will walk you through the process, making as many alterations as you need along the way.

We provide your choice of sterile or non sterile samples to make sure you are totally happy with the finished pack.

Our flexible ordering options (standing orders, call off orders, or make-to-orders) allow stock holding to be tailored to the practice needs.

Why Work With Us?

Quality and Reliability: At Pennine Veterinary, we are committed to providing the highest quality veterinary devices. Our rigorous quality control processes ensure that every product is reliable and safe to use.

Customer Focus: Your satisfaction is our priority. From the moment you place an order to the moment it arrives with you; we strive to provide an exceptional customer experience. Our customer service team is always ready to assist you with any questions or concerns.

Sustainability: We are committed to increasing our sustainable practices in our packaging processes, and our Custom Procedure Packs allow us to do just that. We aim to minimise our environmental impact while providing high-quality products.

Experience: Not only do we have over 60 years of experience within UK manufacturing, but we are also proud to have pioneered the manufacture of Custom Procedure Packs since 1992.

At Pennine Veterinary, we are dedicated to supporting veterinary professionals with high-quality, reliable, and innovative solutions. Our Custom Procedure Packs are all about partnership, and we take responsibility for all aspects of a procedure pack.

Benefits and Specialties

Pennine work in partnership with individual practices/surgeons to create customer specific surgical theatre packs, for a range of specialties;

- ✓ Orthopaedic
- ✓ Soft Tissue
- ✓ Neurology
- ✓ General
- ✓ Cardiology
- ✓ Ophthalmology
- ✓ Laparoscopy

We work closely with veterinarians to design their custom kit's, which can contain a vast range of components such as; swabs, drapes, sutures, blades, tubing, suction devices, wound dressings, bandages, syringes, needles and gowns.

These are then assembled in our cleanroom in Derbyshire, and sterilised locally.

These packs are hugely beneficial for a variety of reasons;

- ✓ Eradication of errors in setting up
- ✓ Reduces emergency surgery response times
- ✓ Reduces lay-up times for procedures
- ✓ Reduces ordering and administration costs
- ✓ Reduces the potential for infection
- ✓ Easier to calculate the on-cost to the customer
- ✓ Ease of use



We're here to help, please get in touch via telephone on +44 (0)1332 794880 or email us at customer.services@penninehealthcare.co.uk



03.

Create Your Bespoke Pack

Here at Pennine, it could not be easier to create your own bespoke pack – in only 4 simple steps. Firstly, let us know which components you require, followed by your drape and gown preferences. We will then produce a sample for you to trial, before final sign-off and order.



We're here to help, please get in touch via telephone on +44 (0)1332 794880 or email us at customer.services@penninehealthcare.co.uk



Standard Surgical Packs

Our pre-made surgical packs are a great option for the practices that don't require large numbers of stock, but want a quick, 'off the shelf' solution to surgical consumables.

These are preassembled, sterile, and readily available to order.

Practices that have not used packs before, may also benefit from these as a good introduction to procedure kits. They can act as a base to customise your own pack.

Orthopaedic Pack

Product Code	Description
PNV-00100	Orthopaedic Pack
1	Trolley cover 152cm x 190cm
1	Bowl 1000ml
1	Connecting tube, soft flexi rib 3m
4	Drape 75cm x 75cm with adhesive edge
1	Extremity drape 100cm x 160cm with 4cm elasticated fenestration
1	Kidney dish 800ml
5	Swabs 7.5cm x 10cm X-RAY
5	Lap swabs 30cm x 30cm X-RAY
1	Bulb Syringe

General Surgery Pack

Product Code	Description
PNV-00400	General Surgery Pack
1	Trolley cover 152cm x 190cm
1	Bowl 1000ml
1	Connecting tube, soft flexi rib 3m
4	Drape 75cm x 75cm with adhesive edge
1	Drape 150cm x 180cm with adhesive edge
1	Kidney dish 800ml
5	Swabs 7.5cm x 10cm X-RAY
5	Lap swabs 30cm x 30cm X-RAY
1	Bulb Syringe



We're here to help, please get in touch via telephone on +44 (0)1332 794880 or email us at customer.services@penninehealthcare.co.uk



Drape Packs

Our standard drape packs are extremely useful as an easy and good practice draping solution for varied animal sizes.

The adhesive edges allow you to create the necessary size of fenestration, without any need to cut into the drapes themselves.

There are four drapes in each pack to achieve the 'four square' drape system. The trolley cover wrapping creates a sterile field so the pack is ready to be used.


These drape packs are great value for money and save on ordering and set up times.

Large Drape Pack

Product Code	Description
VKT-04601	Large Drape Pack
1	Trolley cover 150cm x 190cm, zone reinforced, absorbent
4	Drape 90cm x 90cm with adhesive edge
1	Hand towel 38cm x 40cm

Small Drape Pack

Product Code	Description
VKT-01400	Small Drape Pack
1	Trolley cover, 150cm x 190cm, zone reinforced, absorbent
4	Drape 90cm x 75cm, Bi-laminate with adhesive edge



We're here to help, please get in touch via telephone on +44 (0)1332 794880 or email us at customer.services@penninehealthcare.co.uk

Spay Packs

Our Spay Packs are an excellent solution for spay surgeries as they are great value for money and have all the necessary consumable components in one sterile pouch, ready to use.

One kit allows for less time and waste when opening individual items, and a clear cost.

These are a stocked item, readily available to order.

As on page 3 of this brochure, we also offer fully-customised Packs, allowing for a truly bespoke experience.

Cat Spay Drape Pack

Product Code	Description
VKT-04400	Cat Spay Drape Pack
1	Blue wrap 60cm x 60cm
1	Drape 45cm x 50cm with 3cm x 5cm adhesive fenestration
2	Hand towels 38cm x 40cm
1	Scalpel size 15
5	Swab 10cm x 10cm X-RAY

Bitch Spay Drape Pack

Product Code	Description
VKT-04500	Bitch Spay Drape Pack
1	Blue wrap 60cm x 60cm
1	Drape 100cm x 100cm with 12cm x 9cm adhesive fenestration
2	Hand towels 38cm x 40cm
1	Scalpel size 15
10	Swab 10cm x 10cm X-RAY

Laparoscopic Bitch Spay Pack

Product Code	Description
VKT-07201	Laparoscopic Bitch Spay Kit
1	Drape, patient 300cm x 160cm, bi-laminate, non-reinforced, non-adhesive
1	Trolley cover, 150cm x 190cm, zone reinforced, absorbent
1	Insufflation set, with push fit / male luer connection, 3m x 6mm, square filter
1	Blade, size 11
1	Syringe, 20ml, luer slip
1	Needle, hypodermic, 21g, 1.5 inch, green
1	OR Towel, blue, x-ray, absorbent, 45 x 60cm
1	Bag, breather, 35cm x 50cm, 2fs
1	Suture, biosyn, 3-0, violet, 75cm, 1/2 circle round bodied needle
1	Bag, extras wallet, 18cm x 24cm
10	Swab, gauze, x-ray, 10cm x 10cm, 16ply, without tapes, white, 17 thread





Pennine Veterinary
Ivor Shaw Limited t/a Pennine Healthcare
300 City Gate, London Road, Derby, DE24 8WY, UK.
Tel: +44 (0)1332 794877
www.pennineveterinary.co.uk